Johns Hopkins Engineering

Methods in Neurobiology

Introduction to Behavioral Tests in Neurodegeneration




Examples of behavioral test used in neuroscience and neurodegeneration

Test	Species	Special Considerations		
Basic motor and sensory function				
Homecage activity	Rat, mouse	Can be done acutely or over 24-hr cycle		
Rotarod	Mouse	Requires training		
Hot plate (nociception)	Rat, mouse	Care must be taken to avoid injury		
Learning and memory				
Morris water maze	Rat, mouse	Requires multiple days of training		
Barnes maze	Rat, mouse	Requires multiple days of training		
Radial arm maze	Rat, mouse			
Object recognition	Rat, mouse			
Fear conditioning	Mouse	Can be cued or contextual, requiring different neural circuits; subject to extinction and reinstatement		
Social behavior				
Social interaction/preference	Rat, mouse	Requires a conspecific target, isolation housing is recommended prior to testing		
Sexual behavior	Rat	Should be assessed during dark cycle		
Maternal behavior	Rat, mouse			


Anxiety and depression-like behaviors			
Forced swim test	Rat, mouse	Rats typically require two swim sessions for sufficient immobility	
Tail suspension test	Mouse		
Elevated plus maze	Rat, mouse	Highly sensitive to prior handling/testing	
Light/dark box	Mouse		
	Rat, mouse	Responses vary by age, should be conducted in a soundproof chamber	
Prepulse inhibition	Rat, mouse	Similar to acoustic startle	
Reward			
Conditioned place preference	Rat, mouse	Requires multiple training sessions; subject to extinction and reinstatement	
Sucrose preference	Rat, mouse		
Self-administration	Rat, mouse	Most commonly done in rats; various reinforcers can be used; subject to extinction and reinstatement	
		<u> </u>	

Behavioral test to assay learning and memory in ND

Y Maze





Morris Water Maze



Novel Object Recognition










Training Day

Behavioral test to assay motor function in ND

Open Field task







Cylinder Test

Limb clasping



References

Slide	Reference
2	Teegarden, S. 2020 Behavioral Phenotyping in Rats and Mice. Labome https://www.labome.com/method/Behavioral-Phenotyping-in-Rats-and-Mice.html
3-4	Behavioral and functional neuroscience laboratory. Stanford Medicine. https://med.stanford.edu/sbfnl/services/bm/sm/CylinderTest.html

